




北京爱迪博生物科技 BEIJING ADIPOBIOTECH

Rabbit Anti-Rac1 (Human, Rat, Mouse) Antibody

Catalog No.:	ADB-RAC1-50
Size:	50 µg
Immunogen:	human Rac1 (69-73) (Phe-Leu-Ser-Tyr-Pro)
Host:	rabbit
Antibody Type:	The purified IgG from rabbit serum using an immunogen affinity chromatography.
Formulation:	50 µl of PBS (150mM NaCl, pH 7.4, 0.02% sodium azide and 50% glycerol) in liquid form.
Specificity:	human, mouse and rat Rac1 protein in Western blot.
Applications:	Western blot: 1: 500~ 1000 dilution
Storage:	Store at 4 °C for 2-3 weeks. Store at -20 °C for one year.

订货电话: 010-81786624
010-81786244
Info@AdipoBiotech.com



Detection of Rac1 in extract from A431 cell untreated or treated with EGF (200 ng/ml, 5 min) by Rabbit Anti-Rac1 Antibody